

# UNICHROMAT-Model year 2011

Special Chromatography Cabinet for Äkta-Systems  
from GE-HEALTHCARE

## UNICHROMAT 1500

Mod. 2011 with more power: **1744W**

**New isolation: 80mm**



### UNICHROMAT 1500 for the Äkta System

- Left chamber: Shelf for the Äkta system max 100kg, 2x power socket
- Right chamber: 3x shelves, max. 20kg, 2x power socket
- In the middle: 2x support system for columns.



- Available volumes: 700 / 1500 /1700/2100 liters
- Interior made of stainless steel
- White housing, powder coated
- Digital temperature control with visual alarm
- Digital display in 0,5°C
- Potential-free outlet for connection to a house alarm system,
- **Over-temperature protection**
- **Under temperature safety device with own sensor**
- **The humidity control is based on the evaporator fan functioning according to the compressor working status.**
- CFC/HCFC - free cooling system, fully hermetically sealed and with forced ventilation, mounted on anti-vibration damper, smooth running and low noise emission
- Working temperature: +4°C constant – **with amplified cooling performance, for compensating the heating emission of the Äkta-System**
- Regulation from +1°C to +18°C
- PT-100 temperature control sensor
- Adjustable shelves of stainless steel
- Fully insulated safety-glass doors, self-closing, large glass area for best viewing conditions
- Fully-automatic defroster
- 2x moisture-proofed internal electrical sockets in each chamber
- 2x access ports 40mm left & right
- Illumination
- Double blower for air circulation
- Condensate removal system
- Evaporator is mounted outside of the chamber on top against electromagnetic interference of the Äkta
- Wheels with breaks + 1x adjustable foot=1500/1700/2100
- 4 x adjustable foot=Unichromat 700
- Power 220V / 50Hz
- Lock included

### Specifically designed for the Äkta-System from GE Healthcare:

- Telescopic-shelves, placed according to user's needs
- Retractable over telescopic rails
- Maximum loading capacity 100kg

### Accessories:

- Additional moisture-proofed plugs
- Additional shelves
- Retractable drawers
- Chromatography-column holder
- Data logger with alarm and software
- Models 1500 without middle-column for big instruments
- **Twin Compressor for 2x ÄKTA System (2x1700W)**
- **Special UNICHROMAT for the Äkta AVANT System (Unichromat1700)**

Beste Qualität für das beste Chromatography System.  
Best quality for the best Chromatography System

## UNICHROMAT-1500



## UNICHROMAT-700



**Unichromat 700**  
for the ÄKTA System  
-Shelve for the Äkta System  
-4x Power Socket  
-1x Support System  
for columns

<b>Model</b>	<b>UNICHROMAT 700</b>	<b>UNICHROMAT 1500</b>	<b>UNICHROMAT 2100</b>
<b>Dimensions, mm (WxDxH)</b> <u>external dimensions</u>	<u>780x 840x 2100</u>	<u>1560x 840x 2020</u>	<u>2250x 840x 2020</u>
<b>internal dimensions</b>	600x 670x 1510	<b>1350x 670x 1520</b>	2050x 670x 1520
<b>No. of shelves</b>	3 (without Äkta shelve) 0 (with Äkta shelve)	6 (without Äkta shelve) 3 (with Äkta shelve)	9 (without Äkta shelve) 6 (with Äkta shelve)
<b>WxD, shelve</b>	525x490 mm	525x490 mm	525x490 mm
<b>WxD, Äkta sliding support</b>	550x560 mm	550x560 mm	550x560 mm
<b>No. of glass doors</b>	1	2	3
<b>Weight, net</b>	150 kg	250 kg	350 kg
<b>Temperature range</b>	+1°C bis +18°C	+1°C bis +18°C	+1°C bis +18°C
<b>Power Refrigerant</b>	850W R 404 A (0.35 kg)	<b>1744W</b> R 404 A (0.40 kg)	2407W R 404 A (0.55 kg)

Ihr Fachhändler:

MS-L GmbH \* Kälte- & Temperiertechnik \* Schönbornstr. 7 \* D-69234 Dielheim bei Heidelberg  
Tel: (0)6222 / 93 56-60 \* Fax: (0)6222 / 93 56-66 \* eMail: Info@MS-L.de \* www.MS-L.de