

# D300 Diamond<sup>®</sup> Tips

*High-Quality Performance Tips*



## **Gilson's D300 Diamond Range Means:**

- Pure quality
- Prevention of contamination
- Quality management
- Choice of packaging

**Ideal for molecular biology applications such as PCR, sequencing, and immunology assays.**

Diamond Tips in racks are matched perfectly to Pipetman Ultra Multichannel, thanks to Gilson's patented Rocky Rack™ design.

*Green Signifies ...*

*20 µL to 300 µL*

*... in the Diamond Range*



# D300 Diamond<sup>®</sup> Tips

*Wide Choice of Options and Packaging  
For Different Applications*



## Pure Quality

Made from the purest polypropylene to ensure quality. Certifiable trace metal levels. Tips are untouched by human hands to prevent external contamination (DNA, RNase, etc.).

## Prevention of Contamination

Diamond Tips & Diamond Filter Tips may be autoclaved. Both are also available gamma-ray sterilized, saving you time by eliminating the need to autoclave.

## Quality Management

Color coded for ease of identification—green packaging signifies 20 to 300 µL in the Diamond Range. Fully traceable, with mold number on tips and batch number on bags, boxes, and certificates. Quality certificate in every box.

## Choice of Packaging

- Eco Pack™—The economical solution...box of 10,000 (10 bags of 1000)
- Easy Pack™—Box of 1000 (5 bags of 200)
- Tipack™—Ready-to-use boxes of 96 tips\*
- Tower Pack™—Innovative refill system (10 x 96)\*

\* All of these racked tips benefit from Gilson's Rocky Rack™—dome shaped for comfortable use and perfect fit.



### Description

Eco Pack: D300, box of 10000 (10 bags of 1000 tips)

Easy Pack: D300, box of 1000 (5 bags of 200 tips)

Tipack: Rocky Rack, D300, 10 Tipacks of 96 Diamond Tips in racks

Tipack: D300ST, 10 Tipacks of 96 Sterilized Diamond Tips in racks

Tipack: DF300A, 10 Tipacks of 96 Autoclavable Diamond Filter Tips in racks

Tipack: DF300ST, 10 Tipacks of 96 Autoclavable Diamond Sterilized Filter Tips in racks

Tower Pack: D300, refill system for Diamond Tips

Tower Pack: D300ST, refill system for Sterilized Diamond Tips

### Reference

F161730

F161731

F161732

F161735

F161734

F161733

F167105

F167205

**ISO9001 Certified**

### World Headquarters

#### Gilson, Inc.

3000 W. Beltline Hwy., P.O. Box 620027, Middleton, WI 53562-0027, USA  
Telephone: (1) 800-445-7661 or (1) 608-836-1551 • Fax: (1) 608-831-4451

### Gilson S.A.S.

19, avenue des Entrepreneurs, BP 145, F-95400 VILLIERS LE BEL France  
Telephone: (33-1) 34-29-50-00 • Fax: (33-1) 34-29-50-20

Internet: [www.gilson.com](http://www.gilson.com)

E-mail: [sales@gilson.com](mailto:sales@gilson.com), [service@gilson.com](mailto:service@gilson.com), [training@gilson.com](mailto:training@gilson.com)



LT800439E. Printed in France, December 2002. Specifications subject to change without notice. Gilson, Inc., Gilson S.A.S., and Cyberlab, Inc. are ISO 9001 certified.