

-85/-40°C

**ULTRA LOW
TEMPERATURE
FREEZERS**

ULUF 60/120

ULUF (Ultra Low Undercounter Freezer)

Features on ULUF series

- Single compressor
- Low energy consumption
- Low noise level
- Low heat dissipation
- Insulated inner doors
- Castors and door key lock
- Heated door frame
- Heated vacuum valve
- 100% HCFC/CFC free
- Porthole for external temperature probes


Controller features

- Micro processor controller with digital display
- Approx. 72 hours battery back up for alarms, loggings and temperature display in case of power cut
- Visual and acoustic alarm
- Adjustable high/low temperature alarm
- Power failure alarm
- Probe failure alarm
- Instrument failure alarm
- Open door alarm
- Contact for remote alarm
- Prepared for GSM alarm
- Integrated data logger (software included)
- RS485/232 Interface
- Computer USB data read out
- Auto cycle if probe failure
- Ambient temperature display
- Strows all alarms
- Direct download/upload on/ from USB memory stick
- 3-level password protected
- Battery level indication
- Integrated memory

The ULUF 60/120 is another compact ultra low temperature freezer. With its 55 or 115 Liters capacity, this unit enables you to keep your most used samples right next to you, either in your office or your laboratory room. You will appreciate the easy access and space saving design of this unit, as well as maximum performance and reliability.



ULUF 120



Model	ULUF 60	ULUF 120
Item code (DAI no.)	?	DAI 1417
Dimensions exterior WxDxH (mm)	600x680x810	950x725x810
Net capacity (L)	55	115
Temperature range (°C)	-85/-40	-85/-40
Power supply (V)	230	230
Frequency (Hz)	50/60	50/60

We reserve the right to change specifications without notice. Subject to confirmation and availability.



MS-L GmbH
Kälte- und Temperiertechnik
Schönbornstr. 7
D-69234 Dielheim bei Heidelberg

Tel: (0)6222 / 93 56-60
Fax: (0)6222 / 93 56-66
eMail: Info@MS-L.de